

Mehr Hospital Laboratory

راهنمای نگهدارنده ادرار ۲۴ ساعته

Analyte	یخچال	فریزر	HCL 6mol/L	Boric acid	Acetic acid	Other additives
Albumin	√	√		√		
Aldosterone	√	√	√	√		
Amylase	√					
B2 Microglobuline	√	√				
Calcium			√			
Citrate	√	√	√	√		Toluene
Cortisol	√	√		√	√	
Copper			√			
Creatinine	√	√	√	√		
Cystine		√	√			
Na & K	√			√		
Glucose				√		
5 HIAA	√		√	√		Protection from Sunlight
Immuno electrophoresis	√	√				
17 Ketosteroids			√	√	√	
Mg	√		√			
Metanephrines			√			
Microalbumin	√	√		√		
Normetanephrines		√	√		√	
Osmolality		√				
Oxalate	√		√			
Phosphorous	√	√	√	√		
Protein	√			√		
Urea	√			√		
Uric acid	√	√		√		Toluene or NaHCO3
VMA		√	√	√		
Zink			√			